USPE VET
Veterinary & Animal Care Equipment

Stainless Animal Necropsy Sink Service Station

MODEL USPE-PET-NC-SS | PLANNING & CARE MANUAL

Form USPE-PET-NC-SS-PCM-001 · Rev. A · October 28, 2026

Coordinate an integral sink, continuous rinse, water controls and spray hose in one animal-science service station. Select the eyewash arrangement, drain option and additional receptacle configuration for written review.

Document scope

USPE PET veterinary equipment planning and care reference. The equipment platform is unchanged; the veterinary application requires review. Not a released operating or installation manual.

Plan before you connect

Qualified personnel must verify dimensions, support, service access, utilities, backflow protection and electrical protection. Isolate utilities before maintenance. Do not commission damaged equipment or equipment with leaks, loose supports or electrical defects.

Veterinary application, specimen dimensions and load, component compatibility, drainage, electrical supply and any ventilation requirement must be confirmed in the written quote and approved project documents. This station is not represented as a biosafety cabinet, containment system or clinical fluid-administration device.

USPE PET supplies veterinary necropsy stations, animal-science equipment and coordinated laboratory casework. A division of American Mortuary Coolers.

Equipment and configuration reference

- Integral 304 stainless sink basin, 20□ projection from the cabinet
- Continuous rinse
- Hot and cold gooseneck mixing faucet with paddle / wrist-blade controls
- Hydro aspirator with vacuum breaker and flexible hose
- Spray hose with 8 ft flexible hose
- GFCI duplex receptacle, 120 V
- Removable back cover for plumbing access
- 14-gauge (2 mm) 304 stainless steel, heliarc welded and polished, crevice-free
- Floor or wall mount

Listed specifications

Overall	32″ W × 40″ D × 41″ H
Cabinet depth	20″ (sink basin projects a further 20″)
Sink rim height	23″ from the floor
Clearance under basin	9½″
Overall height	70½″ to the top of the faucet spout
Construction	304 stainless steel, 14 gauge (2 mm), heliarc welded, polished
Faucet	Hot/cold gooseneck mixing faucet, paddle / wrist-blade controls
Aspirator	Hydro aspirator with vacuum breaker and flexible hose; cold water
Spray hose	8 ft flexible hose sprayer; cold water
Electrical	120 V, 60 Hz, 1 phase, 15 A; GFCI duplex receptacle
Water supply	Hot and cold; 58–60 psi and 4 GPM at the hydro aspirator
Rear rough-in	1½″ drain centered 16″ from the left, 7½″ above the floor; hot ½″ at 12″ and cold ½″ at 20″ from the left, both 15″ above the floor; 4″ × 4″ junction box centered 27″ from the left, 21″ above the floor (rough-in sheet)
Mounting	Floor or wall mount

Exact quote configurations

Select the catalog code and listed configuration. A selection is a quotation request, not a purchased item or an installation approval.

USPE PET catalog code	Documented configuration
USPE-PET-NC-SS-01	No Eyewash / Standard Drain / None
USPE-PET-NC-SS-02	No Eyewash / Standard Drain / 1 Additional GFCI
USPE-PET-NC-SS-03	No Eyewash / Standard Drain / 2 Additional GFCI
USPE-PET-NC-SS-04	No Eyewash / 2″ Drain Upgrade / None
USPE-PET-NC-SS-05	No Eyewash / 2″ Drain Upgrade / 1 Additional GFCI
USPE-PET-NC-SS-06	No Eyewash / 2″ Drain Upgrade / 2 Additional GFCI
USPE-PET-NC-SS-07	Eyewash – Left Side / Standard Drain / None
USPE-PET-NC-SS-08	Eyewash – Left Side / Standard Drain / 1 Additional GFCI
USPE-PET-NC-SS-09	Eyewash – Left Side / Standard Drain / 2 Additional GFCI
USPE-PET-NC-SS-10	Eyewash – Left Side / 2″ Drain Upgrade / None
USPE-PET-NC-SS-11	Eyewash – Left Side / 2″ Drain Upgrade / 1 Additional GFCI
USPE-PET-NC-SS-12	Eyewash – Left Side / 2″ Drain Upgrade / 2 Additional GFCI
USPE-PET-NC-SS-13	Eyewash – Right Side / Standard Drain / None
USPE-PET-NC-SS-14	Eyewash – Right Side / Standard Drain / 1 Additional GFCI
USPE-PET-NC-SS-15	Eyewash – Right Side / Standard Drain / 2 Additional GFCI

USPE PET catalog code	Documented configuration
USPE-PET-NC-SS-16	Eyewash – Right Side / 2□ Drain Upgrade / None
USPE-PET-NC-SS-17	Eyewash – Right Side / 2□ Drain Upgrade / 1 Additional GFCI
USPE-PET-NC-SS-18	Eyewash – Right Side / 2□ Drain Upgrade / 2 Additional GFCI
USPE-PET-NC-SS-19	Eyewash + Drench Hose – Left Side / Standard Drain / None
USPE-PET-NC-SS-20	Eyewash + Drench Hose – Left Side / Standard Drain / 1 Additional GFCI
USPE-PET-NC-SS-21	Eyewash + Drench Hose – Left Side / Standard Drain / 2 Additional GFCI
USPE-PET-NC-SS-22	Eyewash + Drench Hose – Left Side / 2□ Drain Upgrade / None
USPE-PET-NC-SS-23	Eyewash + Drench Hose – Left Side / 2□ Drain Upgrade / 1 Additional GFCI
USPE-PET-NC-SS-24	Eyewash + Drench Hose – Left Side / 2□ Drain Upgrade / 2 Additional GFCI
USPE-PET-NC-SS-25	Eyewash + Drench Hose – Right Side / Standard Drain / None
USPE-PET-NC-SS-26	Eyewash + Drench Hose – Right Side / Standard Drain / 1 Additional GFCI
USPE-PET-NC-SS-27	Eyewash + Drench Hose – Right Side / Standard Drain / 2 Additional GFCI
USPE-PET-NC-SS-28	Eyewash + Drench Hose – Right Side / 2□ Drain Upgrade / None
USPE-PET-NC-SS-29	Eyewash + Drench Hose – Right Side / 2□ Drain Upgrade / 1 Additional GFCI
USPE-PET-NC-SS-30	Eyewash + Drench Hose – Right Side / 2□ Drain Upgrade / 2 Additional GFCI

Casework and room handoff

Schedule the station with its base, wall and tall cabinets. Mark quantity, room zone, cutouts, handedness, height, support, projections and service clearances. Confirm finishes, glass and hardware in the quote.

Casework material, glass and door-style form

Stainless steel laboratory casework is our specialty. Choose standard stainless steel, powder-coated stainless steel or powder-coated steel. The 100% American-made casework message does not establish the origin of every purchased fixture or accessory.

Cabinet / room zone / quantity: _____
Material / grade / gauge: _____
Body color & tier / front color & tier: _____
Glass appearance: clear / frosted / grey-tinted / bronze-tinted / not applicable
Door style: solid hinged / glazed hinged / sliding / open storage
Handing / frame / glass thickness / safety-glazing requirement: _____
Handle / hinge / sample / drawing revision: _____

Glass appearance is a request, not a stocked glazing specification or separate SKU. Confirm tempered or laminated requirements, thickness, edge treatment, hardware and cabinet fit in the written quote. No impact, fire, containment or cleanroom rating is implied.

Receiving and commissioning

Check condition, model and options against the approved schedule. Report damage or missing components before use. Obtain component instructions, complete utility checks and document commissioning. A photograph does not establish capacity, containment or waste-disposal approval.

Care and service

Follow component instructions and the facility's approved cleaning protocol. Confirm agent, dilution, contact time and material compatibility. No pathogen-inactivation or sterilization claim is made. Do not mix chemicals or clean through electrical components.

Stop use if damage, leaks or electrical defects occur; qualified personnel must isolate and assess them. Order replacements by exact model and component details, not appearance.

Documentation and quote support

1-888-792-9315 · po@mymortuarycooler.com. Request current model-specific operating instructions, approved drawings and applicable warranty terms with your equipment schedule.

https://www.mymortuarycooler.com/pages/uspe-quote-request?brand=uspe-vet