USPE VET
Veterinary & Animal Care Equipment

Wall-Mount Animal Necropsy Service Station

MODEL USPE-PET-NC-WM | PLANNING & CARE MANUAL

Form USPE-PET-NC-WM-PCM-001 · Rev. A · October 28, 2026

Coordinate a wall-mounted utility and sink position for a veterinary necropsy room. Keep the supporting wall, service access and fixture arrangement on the approved room drawing.

Document scope

USPE PET veterinary equipment planning and care reference. The equipment platform is unchanged; the veterinary application requires review. Not a released operating or installation manual.

Plan before you connect

Qualified personnel must verify dimensions, support, service access, utilities, backflow protection and electrical protection. Isolate utilities before maintenance. Do not commission damaged equipment or equipment with leaks, loose supports or electrical defects.

Veterinary application, specimen dimensions and load, component compatibility, drainage, electrical supply and any ventilation requirement must be confirmed in the written quote and approved project documents. This station is not represented as a biosafety cabinet, containment system or clinical fluid-administration device.

USPE PET supplies veterinary necropsy stations, animal-science equipment and coordinated laboratory casework. A division of American Mortuary Coolers.

Equipment and configuration reference

- Heavy-duty 304 stainless steel for corrosion resistance
- Spray-hose assembly with cold-water control, nozzle, and 10-ft flexible hose
- 0.75 HP waste disposer for efficient waste removal
- Deluxe reverse-flow aspirator that helps prevent suction clogs
- Hot and cold water unit with lever-handle operation
- Utility faucet with swivel spout
- Built-in GFCI duplex receptacle
- Open, accessible back for plumbing access and modifications
- Optional emergency eyewash attachment

Listed specifications

Height	71"
Length (left to right)	30"
Depth (front to back)	24.5"
Sink basin	12" × 16" × 8"
Material	304 stainless steel
Power	110V with GFCI-protected receptacle

Exact quote configurations

Select the catalog code and listed configuration. A selection is a quotation request, not a purchased item or an installation approval.

USPE PET catalog code	Documented configuration
USPE-PET-NC-WM-01	Standard Drain
USPE-PET-NC-WM-02	2□ Drain Upgrade

Casework and room handoff

Schedule the station with its base, wall and tall cabinets. Mark quantity, room zone, cutouts, handedness, height, support, projections and service clearances. Confirm finishes, glass and hardware in the quote.

Casework material, glass and door-style form

Stainless steel laboratory casework is our specialty. Choose standard stainless steel, powder-coated stainless steel or powder-coated steel. The 100% American-made casework message does not establish the origin of every purchased fixture or accessory.

Cabinet / room zone / quantity: _____
Material / grade / gauge: _____
Body color & tier / front color & tier: _____
Glass appearance: clear / frosted / grey-tinted / bronze-tinted / not applicable
Door style: solid hinged / glazed hinged / sliding / open storage
Handing / frame / glass thickness / safety-glazing requirement: _____
Handle / hinge / sample / drawing revision: _____

Glass appearance is a request, not a stocked glazing specification or separate SKU. Confirm tempered or laminated requirements, thickness, edge treatment, hardware and cabinet fit in the written quote. No impact, fire, containment or cleanroom rating is implied.

Receiving and commissioning

Check condition, model and options against the approved schedule. Report damage or missing components before use. Obtain component instructions, complete utility checks and document commissioning. A photograph does not establish capacity, containment or waste-disposal approval.

Care and service

Follow component instructions and the facility's approved cleaning protocol. Confirm agent, dilution, contact time and material compatibility. No pathogen-inactivation or sterilization claim is made. Do not mix chemicals or clean through electrical components.

Stop use if damage, leaks or electrical defects occur; qualified personnel must isolate and assess them. Order replacements by exact model and component details, not appearance.

Documentation and quote support

1-888-792-9315 · po@mymortuarycooler.com. Request current model-specific operating instructions, approved drawings and applicable warranty terms with your equipment schedule.

https://www.mymortuarycooler.com/pages/uspe-quote-request?brand=uspe-vet